ADDRE

Supporting Health

CATALOG 2023



Pain Management Aids
Rehabilitation Aids
Immobilizers
Traction Kits

www.AdoreMedical.com



Our family stands behind Adore Medical so you can stand up for yours.

-**Damodar Prasad Agrawal**Founder



Our mission is to make high quality orthopaedic support braces, fracture aids and pain management aids easily accessible because we believe that everyone deserves good health.

We strictly follow international standards and source only the highest quality materials to develop our products.

Our commitment to supporting good health for everyone drives every facet of our business and partnerships. Innovation, dedication and strong relationships with our distribution partners are the core values that define Adore Medical Private Limited.



CONTENTS

05

ABDOMEN & SPINE SUPPORT



Abdominal Belt

Tummy Trimmer

Maternity Belt

Ash Brace (Hyper Extension Brace)

Rib Belt

Sternal Support

Lumbo Sacral Support

Contoured L.S. Support

Dorso Lumbar Brace

Hernia Belt

Waist Belt (Neoprene)

Scrotal Support

Posture Corrector

Back Rest

Coccyx Cushion Seat

09

CERVICAL SUPPORT



Cervical Collar Hard

Cervical Collar Soft

Cervical Collar Soft with Support

Cervical Pillow (Memory Foam)

Cervical Pillow Regular

Cervical Traction Kit + Weight Bag

(Sleeping)

Cervical Orthosis

(Philadelphia Collar)

11

SHOULDER & CLAVICLE SUPPORT



Shoulder Support (Neoprene)

Elastic Shoulder Immobilizer

Universal shoulder Immobilizer

Clavicle Brace

Clavicle Brace With Buckle

13

ARM & ELBOW SUPPORT



Pouch Arm Sling (Tropical)

Pouch Arm Sling (Baggy)

Elbow Support

Tennis Elbow Support

Forearm Splint

Tourniquet

15

WRIST, FINGERS & THUMB SUPPORT



Wrist Brace (Neoprene)

Wrist Brace

Wrist Brace with Thumb Support

Wrist Brace With Thumb Support (Neoprene)

Wrist & Forearm Splint (Right/Left)

Elastic Wrist Splint (Right/Left)

Finger Cot

Finger Extension Splint

Frog Splint

Mallet Finger Splint

Thumb Spica Splint

Exercise Gel Ball (Soft/Medium/Hard)

Weight Cuff

19

PELVIC & HIP SUPPORT



Pelvic Binder

Pelvic Traction Kit + Weight Bag

Traction Weight Bag

Leg Traction Brace

20

THIGH, KNEE, ANKLE & CALF SUPPORT



Skin Traction Set (PUF Liner)

Thigh Support

Ankle Brace

Ankle Binder

Anklet

Knee Cap

Knee Cap Open Patella

Knee Cap Rigid Hinge

Knee Wrap (Neoprene)

Knee Wrap Hinged (Neoprene)

Elastic Knee Support

Patellar Support

Compression Stocking Mid Thigh

Compression Stocking Below Knee

Medical Compression Stocking

Mid Thigh - Classic

Medical Compression Stocking

Below Knee - Classic

Knee Immobilizer (14")

Knee Immobilizer (19")

Knee Immobilizer (22")

R.O.M. Knee Brace (22")

Calf Support

26

FOOT, HEEL & TOES SUPPORT



Foot Drop Splint (Right/Left)

Arch Support

Toe Separator

Insole (silicone)

Heel Cushion (silicone)

Metatarsal Pad

Cast Shoe

28

HOT & COLD THERAPY AIDS



Hot & Cold Pack

29

PHYSIO AIDS



Exercise Band (Actiband)

30

OTHER PRODUCTS

Disposable 3 Ply Face Mask

Cotton Crepe Bandage B.P. (Acticrepe)



ABDOMEN & SPINE SUPPORT





AD-101

Abdominal Belt

Our abdominal belt is perfect for postoperative conditions of the abdomen and post-delivery umbilical and ventral hernias.

Indications:

- Postoperative conditions of the abdomen.
- Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the navel



ΔΠ-10

Tummy Trimmer

Optimal compression of our Tummy Trimmer binds surgical incisions and helps heal abdomen and chest.

Indications:

- Postural deformity.
- Osteoporosis of the spine.
- Postoperative care of the spine.
- Mild thoraces and lumbar injuries.
- Degenerative Disc Disease
- Herniated Disc or PIVD.

Size	Inches	Cm
S	28-32	70-80
М	32-36	80-90
Ĺ	36-40	90-100
XL	40-44	100-110
XXL	44-48	110-120
XXXL	48-52	120-130
XXXXL	52-56	130-140



*Measure circumference around the navel



AD-111

Maternity Belt

This maternity belt provides relief of hip and pelvic pain through sufficient stretching and gentle compression to reduce bladder pressure and spread the baby's weight.

- Postoperative conditions of the abdomen.
- Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
S	32-38	80-95
М	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3' below the navel





AD-108

Ash Brace (Hyper Extension Brace)

Our Ash Brace is aerodynamically shaped with telescoping, aluminum section, and malleable horizontal main body to form a cross bar.

Siz	ze	Inches	Cm	
	Short	Up to 5.6 Ft	Up to 165	H
Universal	Long	Above 5.6 Ft	Above 165	

*Measure height of the patient

Indications:

- Osteoporotic Bone Collapse.
- Ankylosing Spondylitis.
- Conservative management of Vertebral Bone Tuberculosis.
- Postoperative, post-traumatic, or chronic spinal pathology.
- Complications of thoracolumbar junction which needs moderate Immobilization.



Rib Belt

Designed with a splinting pad for comfort, our Rib Belt is useful in aiding the healing process in rib fracture, relieving chest pain, and intercostal muscle strain.

Indications:

- Postoperative conditions of the abdomen.
- · Post-delivery.
- Umbilical and ventral hernia.

Size	Inches	Cm
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145

*Measure circumference around the chest below underarms



Sternal Support

Designed for comfort with an easy hook & loop closure, moisturewicking fabric and splinting pad.

Indications:

- · Post cardiothoracic surgeries.
- Fracture of the sternum.
- Pain due to chest injury.
- · Rib fractures.
- Osteoporosis.

Size	Inches	Cm
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145

*Measure circumference around the chest below underarms



AD-102

Lumbo Sacral Support


Our lumbo-sacral support is used to comfortably maintain correct posture and support the lower back while sitting, walking or driving.

- Postoperative conditions of the abdomen.
- · Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
Size	Inches	CIII
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145

*Measure circumference around the navel





AD-103

Contoured L.S. Support

This Contoured Lumbo-Sacral Support is ideal for postoperative care of lumbar spine and for lower back pain, strain or sprains.

Indications:

- Postoperative care of the lumbar spine.
- · Low back pain.
- Low back sprains/strains.
- · Lumbar disc injury.
- Degenerative Disc Disease.
- Herniated Disc or PIVD.

Size	Inches	Cm
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the navel



Dorso Lumbar Brace

Features anatomically contoured splints for precise immobilization and clavicle straps with cotton stockinet to correct shoulder posture alignment.

Indications:

- Postural deformity.
- Osteoporosis of the spine.
- Postoperative spine care.
- Mild lumbar and thoraces injury.
- Degenerative Disc Disease.
- Herniated Disc or PIVD.

Size	Feet	Cm
Short	4.6-5.6 Ft.	140-171
Long	5.6-6.6 Ft.	171-201



*Measure height of the patient

Size	Inches	Cm
Universal	28-44	71-112
Special	44-52	112-132



*Measure circumference around the navel



AD-109

Hernia Belt

Designed for a comfortable and supportive fit our hernia belt is lightweight and features adjustable pelvic and leg straps.

Indications:

- Reducible inguinal or scrotal hernia.
- Post hernia repair surgery.
- Post-delivery umbilical and ventral hernia.
- Post abdominal hernia repair surgery.

Size	Inches	Cm	
S	28-34	70-85	
М	34-40	85-100	
Ĺ	40-46	100-115	
XL	46-52	115-130	



*Measure circumference 3' below navel



AD-112

Waist Belt (Neoprene)

Our wrap-around design waist belt is made with neoprene bonded fabric for four-way stretchability and maximum support.

- Weak abdominal muscles.
- Swelling & stiffness in abdominal muscles.
- Postoperative, post-delivery support.

Size	Inches	Cm	
Universal	28-50	70-125	

*Measure circumference around the navel





AD-1103

Scrotal Support

Fabricated extra-soft and hypoallergenic materials, our scrotal support is well-ventilated to allow for ease of use and extended wear comfort.

Indications:

- · Palliative care.
- Treatment for postoperative care.
- Tenderness, pain or swelling in scrotal area.

Size	Inches	Cm	
S	28-34	70-85	1_
М	34-40	85-100	4-
L	40-46	100-115	
XL	46-52	115-130	

*Measure circumference 3" below navel



AD-110

Posture Corrector

Designed to correct slouching and improve posture.

Indications:

- Back pain.
- Spine misalignment.
- Lumbar pain.
- Poor posture.

Size	Inches	Cm
СН	22-28	56-71
S	28-34	71-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130

*Measure circumference around the chest below underarms



AD-110

Back Rest

Promotes good posture and correct spinal alignment by allowing the spine to follow natural curve.

Indications:

- Postoperative discomfort.
- Post-injury pain.
- Poor posture.
- Lower back pain.
- Pregnancy discomfort.

Size	Cm
Universal	152-192



AD-1102

Coccyx Cushion Seat

The high-density foam of our coccyx cushion features a U-shaped cut that gently suspends your coccyx to relieve pressure on your spinal cord.

- Pudendal nerve entrapment pain.
- Chronic pelvic pain syndrome.
- Rectal surgery.
- Prostatitis.
- Anal surgery.

Size	Inches (LXBXH)	Cm (LXBXH)
Universal	16.4X16.5X2.7	41.7X42X7

- Coccyx fracture.
- · Hemorrhoids.
- Interstitial Cystitis.
- Pelvic floor pain.
- Pressure sores.



CERVICAL SUPPORT AIDS





Cervical Collar Hard

Our hard Cervical Collar is made of high-quality polyethylene with padded edges moulded with Ethafoam.

Indications:

- · Cervical Spondylitis.
- Cervical Disc Disease.
- Fracture of cervical spine.
- · Whiplash injuries.

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
М	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25



*Measure distance between edge of sternum & lower jaw with neck in the extension position



Cervical Collar Soft

Made of soft density foam for gentle and cushioned neck support.

Indications:

- Cervical Trauma.
- Neck Pain and Stiffness.
- Excessive Mobility.
- Rheumatoid Arthritis.
- Strain of Neck Muscles.
- · Whiplash Recovery.
- Degenerative Disc Disease.

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
М	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50



*Measure distance between edge of sternum & lower jaw with neck in the extension position




AD-201

Cervical Collar Soft With Support

Made of firm density foam covered with a soft, smooth cotton stockinet splint inside for better neck support.

- Cervical Trauma.
- Neck Pain and stiffness.
- Excessive mobility.
- Rheumatoid Arthritis.
- Strain of neck muscles.
- Whiplash recovery.
- Degenerative disc disease.

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
М	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50



*Measure distance between edge of sternum & lower jaw with neck in the extension position





AD-203

Cervical Pillow (Memory Foam)

Made of memory foam to provide proper support to the neck region while relaxing or sleeping.


Indications:

- Temporomandibular disorders.
- Cervical trauma.
- · Morning hand stiffness and swelling.
- Whiplash injuries.
- Arthritis.
- · Neck pain.

Size	Inches	Cm
Universal	3.0-4.5	7.50-11.25



*Measure distance between edge of sternum & lower jaw with neck in the extension position



Cervical Pillow Regular

Provides sustained support to cervical vertebrae during sleep.

Indications:

- Temporomandibular disorders.
- · Cervical trauma.
- Morning hand stiffness and swelling.
- Whiplash injuries.
- · Arthritis.
- Neck pain.

Size	Inches	Cm
Universal	3.0-4.5	7.50-11.25



*Measure distance between edge of sternum & lower jaw with neck in the extension position



AD-206

Cervical Traction Kit + Weight Bag (Sleeping)

With a width adaptation from 25mm to 50mm, our Cervical Traction Kit is developed to fit on the top of bed.

Indications:

- · Cervical Spondylitis.
- Cervical Disc Disease.
- Fractures of cervical spine.
- Mechanical stress due to poor posture.

Universal



AD-204

Cervical Orthosis (Philadelphia Collar)

Made of Ethafoam material to padded support our Philadelphia Collar has a trachea opening in front.

Indications:

- · Cervical trauma.
- · Cervical sprain.
- · Cervical Spondylitis.

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
М	3.5-4.0	8.75-10
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50




*Measure distance between edge of sternum & lower jaw with neck in the extension position



SHOULDER & CLAVICLE SUPPORT





Shoulder Support (Neoprene)

Made from Neoprene bonded nylon fabric which is four-way stretchable with an anatomical shape to ensure perfect fitting and optimal immobilization.

-	0	N 1 0	at	-	2
-		111	-		115
ш		$a \circ \circ$	u	\cdot	

- Shoulder Injury.
- Rotator cuff injury.
- Shoulder dislocation
- Shoulder inflammation.
- · Arthritic pain of the shoulder joint.
- · Bursitis and Tendonitis.

Elastic Shoulder

Size	Inches	Cm
Universal	26-40	65-100
Special	40-52	100-130

*Measure circumference around the chest below underarms



Immobilizer

Our Elastic Shoulder Immobilizer is made from breathable and moisture transporting weaving fabric.

Indications:

- Shoulder dislocation.
- Post-surgical immobilization.
- Shoulder joint Pain.
- Forearm pain.
- Shoulder sprain.

Size	Inches	Cm
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the chest below underarms




AD-303

Universal shoulder Immobilizer

Our Universal Shoulder Immobilizer comes in a flexible size to support shoulder immobilization.

- Shoulder dislocation.
- Post-surgical immobilization.
- Shoulder joint Pain.
- · Forearm pain.
- Shoulder sprain.





*Measure circumference around the chest below underarms





AD-304

Clavicle Brace

Eight-shaped band design of Clavicle Brace is cushioned from soft polyurethane foam padding.

Indications:

- Clavicle fracture.
- Clavicle strain.
- Clavicle trauma.
- Poor posture problems.

Size	Inches	Cm
СН	22-28	56-71
S	28-34	70-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	132-147


*Measure circumference around the chest below underarms



Clavicle Brace with Buckle

Eight-shaped band design of Clavicle Brace is cushioned from soft polyurethane foam padding.

- Clavicle fracture.
- Clavicle strain.
- Clavicle trauma.
- Poor posture problems.



*Measure circumference around the chest below underarms



ARM & ELBOW SUPPORT





AD-301

Pouch Arm Sling (Tropical)

Our Pouch Arm Sling provides comfortable support to the wrist and arm with an adjustable strap with a buckle for proper immobilization.

Indications:

- Arm fracture.
- · Wrist fracture.
- Arm & Wrist sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast.

Size	Inches	Cm
СН	22-28	56-71
S	28-34	71-85
М	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-147

*Measure circumference around the chest below underarms



AD-306

Pouch Arm Sling (Baggy)

Our baggy Pouch Arm Sling provides light-weight support to the wrist and arm through a thumb cradle, and an adjustable strap with a buckle for proper immobilization.

Indications:

- Arm Fracture.
- Wrist Fracture.
- Arm & Wrist Sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast rehabilitation.

Size	Inches	Cm	
СН	22-28	56-71	
S	28-34	71-85	
М	34-40	85-100	-
L	40-46	100-115	
XL	46-52	115-130	· '
XXL	52-58	130-147	

*Measure circumference around the chest below underarms



AD-504

Elbow Support

Our Elbow Support is made from a soft, elastic material. It provides compression and support for injured tendons.

- Forearm/Elbow pain.
- Post Elbow
- Arthroscopy.
- Sports and workout stress.
- Post-cast rehabilitation.
- Postoperative weakness.
- Arthritis / Bursitis.

Size	Inches	Cm	
S	7-10	18-25	
М	10-12	25-30	Sold Sold Sold Sold Sold Sold Sold Sold
L	12-14	30-35	
XL	14-16	35-40	

*Measure circumference around the elbow in extended position





۸D-E0Z

Tennis Elbow Support

Our Tennis Elbow support has adjustable straps and a reverse fastener for a strong grip and custom compression.

Indications:

- Tennis/Golf Elbow.
- Tendonitis.
- Swelling.

Size	Inches	Cm
S	7-10	18-25
М	10-12	25-30
L	12-14	30-35
XL	14-16	35-40

*Measure circumference around the elbow in extended position



Forearm Splint

Our Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb.

Indications:

- Post-cast/Post-surgery protection and rehabilitation.
- Conservative management of Kienbock's disease.
- Rheumatic disease of the wrist.
- · Arthritic hand positioning.

Size	Inches	Cm
Universal	7.5-10.50	19-26
Special	10.50-13.50	26-34



*Measure distance from wrist to elbow



AD-1201

Tourniquet (Upper & Lower Extremities)

Our Tourniquet is padded with soft foam on the inner surface and constructed with a lasting hook & loop strap with a tough polypropylene fastener.

Size	Inches	Cm	
Universal	9.8-11.10	24.90-28.20	2

*Measure circumference around the elbow in extended position

- Stop constant bleeding.
- Restrain the circulation of blood to an extremity for a period of time.

WRIST FINGERS & THUMB SUPPORT





AD-508

Wrist Brace (Neoprene)

Created out of triple-layered neoprene. The top layer features continuous loops with rugged nylon fabric in the inner layer.

Indications:

- Tendonitis.
- Tenosynovitis.
- Arthritic care of wrist.
- Sprains and strains.
- Dislocation of the wrist joint.
- Postoperative and postcast care.
- Prophylaxis for exercise and sports.

Size	Inches	Cm
Universal	5.0-9.5	12.5-23.80



*Measure circumference evenly around wrist joint



AD-501

Wrist Brace

Designed with adjustable straps, moisture-wicking fabric and durable nylon reinforcement.

Indications:

- Sprains and strains.
- Mild Osteoarthritis.
- Swelling & inflammation.
- Prophylaxis for exercise and sports.

Size	Inches	Cm
S	4.9-6.1	12.5-15.5
М	6.1-7.3	15.5-18.5
L	7.3-8.5	18.5-21.5



*Measure circumference evenly around wrist joint



AD-502

Wrist Brace with Thumb Support

Designed with adjustable wrap-around straps, moisture-wicking fabric and durable nylon reinforcement.

- Sprains and strains .
- Mild Osteoarthritis.
- Swelling & inflammation.
- Tendonitis.

Size	Inches	Cm	P
Size	Inches	CIII	
Universal	5.0-9.5	12.5-23.80	

*Measure circumference evenly around wrist joint





AD 500

Wrist Brace with Thumb Support (Neoprene)

Neoprene based design with adjustable wraparound straps, moisture-wicking fabric and durable nylon reinforcement.

Indications:

- Ligamentous laxity of the thumb.
- Median nerve injury.
- Skier's thumb.
- Carpometacarpal joint injury.
- Post operative and post-cast care.
- · Prophylaxis for exercise and sports.

Size	Inches	Cm
Universal	5.0-9.5	12.5-23.80



*Measure circumference evenly around wrist joint



Wrist and Forearm Splint (Right/Left)

Our Wrist & Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb.

Indications:

- Post-cast management.
- Post-surgical protection and rehabilitation.
- Conservative management of Kienbock's disease.
- Rheumatic disease of the wrist.
- Arthritic hand positioning.

Size	Inches	Cm
S	7.2-8.5	18.3-21.6
М	8.5-9.8	21.6-24.9
L	9.8-11.1	24.9-28.2
XL	11.1-12.4	28.2-31.5



*Measure the circumference at a distance 6" from the wrist along the arm



AD-506

Elastic Wrist Splint (Right/Left)

Our athletic-inspired elastic splint is a lightweight wrist support with double straps and anatomical design for rigid immobilization, while maintaining the wrist in a functional position.

- · Carpal Tunnel Syndrome.
- Post-surgery.
- Post-cast management.
- · Wrist sprains and strains.
- Tendonitis.

Size	Inches	Cm
S	4.9-6.1	12.5-15.5
М	6.1-7.3	15.5-18.5
L	7.3-8.5	18.5-21.5
XL	8.5-9.7	21.5-24.5



*Measure the circumference at a distance 6" from the wrist along the arm





AD-601

Finger Cot

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort.

Indications:

- Finger injury.
- Post finger Fracture.
- Postoperative care.
- · Sprains and bruises of fingers.

Size	Inches	Cm
S	< 2	< 5
М	2-3	5-7.5
L	3-4	7.5-10



*Length of fingers

|--|

AD-604

Finger Extension Splint

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort.

Indications:

- Metacarpal injury.
- Fracture.
- Post-surgery.
- Sprains and bruises of fingers.

Size	Inches	Cm
S	< 2	< 5
М	2-3	5-7.5
L	3-4	7.5-10



*Length of fingers



AD-602

Frog Splint

Made of light-weight but rigid aluminum frame to rehabilitate and align the post ligamentary damage of the Interphalangeal joints.

Indications:

- Interphalangeal (IP) joint Injury.
- Finger Injury.
- Finger sprain.
- Postoperative care.

Size	Inches	Cm
S	< 2	< 5
М	2-3	5-7.5
L	3-4	7.5-10




*Length of fingers

Mallet Finger Splint

Ethafoam liner for cushioning and also suitable for post-operative care.

- Bruising of the dorsum of the DIPJ.
- Mallet finger deformity.
- May help manage the extensor tendon injury.
- Minor fractures of the distal phalanges.
- Distal phalange postfracture care.
- Distal interphalangeal (DIP) joint conditions.
- Swollen fingertip.
- · Post-surgical care.





*Length of fingers



AD-005





AD-603

Thumb Spica Splint

A flexible and comfortable splint made of malleable material for easy application and removal.

Indications:

- Bursitis.
- Sprains of thumb.
- Avulsion fractures of the thumb.
- · Soft tissue injuries.
- Thumb pain symptoms.
- CMC & MCP joint conditions.

Size	Inches	Cm	SA
Universal	5-9.5	12.5-23.8	

*Measure circumference evenly around wrist joint



Exercise Gel Ball

Anti-stress exercise gel ball in three colours depending on hardness. Good for building strength of fingers, palm and wrist.

Indications:

- Post fracture.
- Stress.
- · Weak forearm and grip strength.

Soft	Medium	Hard



AD-901

Weight Cuff

Our weight cuff is fit for graded exercise therapy (GET) and also cosmetic arm muscle shaping.

- Graded exercise therapy (GET).
- Poor muscle tone.
- Post-cast care.

Size	½ Kg	1 Kg	2 Kg
Variations of Weight Cuff fit all wrists and ankles.			

PELVIC & HIP SUPPORT






Pelvic Binder

Our Pelvic Binder is made with ruggedness for extended and repetitive uses.

Indications:

- Lumbar spondylitis.
- · Lower back pain.
- Postural fatigue.
- · Postural deformities.
- · Vertebral fractures.
- · Mechanical stress due to poor posture.

Size	Inches	Cm
S	32-38	80-95
М	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3" below navel



Pelvic Traction Kit + Weight Bag

A versatile traction kit that can be used in all kinds of tractions and on all kinds of beds.

Indications:

- Lumbar spondylitis.
- Lower back pain.
- Postural fatigue.
- Postural deformities.
- Vertebral fractures.
- · Mechanical stress due to poor posture.

Size	Inches	Cm
S	32-38	80-90
М	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3" below navel



AD-803W

Traction Weight Bag

The Weight Bag applies uniform and controlled traction to the lumbar, sacral, and thoracic areas of the vertebral column.

Indications:

- Convenient home use
- Controlled continuous traction
- Easy removal and application


Weight Capacity		
Liquid	10 Ltr.	
Solid	12 Kg.	



Four sets of nylon based buckle closure with an additional locking mechanism to ensure a secure fit.

- Hip/Sacral pain.
- Shaft fractures of the femur.
- Traction to relieve pain in the lower back
- When pelvic traction is not possible or not desired.

Size	Inches	Cm
S	8.5-11.0	21.50-27.90
М	11.0-13.50	27.90-34.20
L	13.5-16.0	34.20-40.60
XL	16.0-18.50	40.60-46.90



*Measure circumference around calf approx 6" below knee



AD-804



THIGH, KNEE, ANKLE & CALF SUPPORT





AD-805

Skin Traction Set (Puf Liner)

Our Skin Traction kit is made with foam liner designed to ensure there is no friction damage to the skin and no slippage during traction.

Indications:

- Need for indirect traction to extremities.
- · Pains & spasms.

Size	Inches	Cm
Child	< 12	< 30
Adult	> 12	> 30

FDORE


Thigh Support

The weaving method for our Thigh Support provides seamless stretching stockings with a tapered top to bottom design to hold in a natural shape.

Indications:

- Hamstring strain.
- · Muscle spasms or cramping.
- Arthritis.
- Stiffness or trauma.
- Pain & inflammation.
- Sports Injury.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6-22	49-55
XL	22-24.4	55-61



*Measure circumference around mid thigh approx 6" above knee



AD-412


AD-416

Ankle Brace

Easy to wear and remove brace made of rigid polypropylene material for strong immobilization.

- Mild ankle sprains
- Recurring ankle injury.
- Tarsal tunnel syndrome.
- Swollen ankle. (Not indicated for severe ankle instabilities)

Size	Inches	Cm
СН	Below 7.2	Below 18
S	7.2-8.4	18-21
М	8.4-9.6	21-24
L	9.6-10.8	24-27
XL	10.8-12	27-30



*Measure circumference approx 2" above the ankle joint





AD-401

Ankle Binder

The figure of eight design makes it comfortably wearable, even with shoes.

Indications:

- Supination injuries.
- Post-supination injuries.
- · Syndesmotic ligament injury.
- Chronic ligament instability in the upper & lower ankle.

Size	Inches	Cm	
S	7.2-8.4	18-21	
М	8.4-9.6	21-24	9-
L	9.6-10.8	24-27	É
XL	10.8-12	27-30	

*Measure circumference approx 2" above the ankle joint



Anklet

Provides compression for the foot and ankle joint.

Indications:

- · Lymphoedema.
- · Lipedema.
- · Chronic venous insufficiency.
- Varicose veins.
- Post Thrombotic Syndrome.
- · Venous Stasis Ulcer.
- Post Sclerotherapy.
- · Lipodermatosclerosis.

Size	Inches	Cm
S	7.2-8.4	18-21
М	8.4-9.6	21-24
L	9.6-10.8	24-27
XL	10.8-12	27-30

*Measure circumference approx 2" above the ankle joint



AD-403

Knee Cap

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Knee Cap Open Patella

Spandex-based fabric ensure

effective compression with

stretchability and comfort.

- · Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees


Size	Inches	Cm	
S	14.8-17.2	37-43	1
М	17.2-19.6	43-49	1
L	19.6-22.0	49-55	6"
XL	22.0-24.4	55-61	1 + 10
XXL	24.4-26.8	61-67	
XXXL	26.8-29.2	67-73	

*Measure circumference around mid thigh approx 6" above knee



Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees



*Measure circumference around mid thigh approx 6" above knee





AD-405


Knee Cap Rigid Hinge

An atomically tapered form and spandexbased fabric ensure effective compression with stretchability and comfort.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- · Quad cramps.
- Knee Sprain.
- · Weak knees.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee

ROORE

AD-415

Knee Wrap (Neoprene)

Made with triple-layer neoprene fabric. The triple hook loop closing system ensures easy wearing, removal and proper fit.

Indications:

- · Arthritic pain.
- Osteoarthritis of the knee.
- · Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees.

Size	Inches	Cm
Universal	14.8-22	37-55
Special	22-29.2	55-73

*Measure circumference around mid thigh approx 6" above knee



AD-410

Knee Wrap Hinged (Neoprene)

Simple and quick put on with an easy wrap-around application made from tapered, four-way stretchable neoprene material.

Indications:


- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.

Elastic Knee Support

Made from elastic, comfortable and durable material with adjustable strap for proper compression and

- · Quad cramps.
- Knee sprain.
- · Weak knees.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67




*Measure circumference around mid thigh approx 6" above knee



customized fit.

- Indications:Arthritic pain.
- Osteoarthritis of the knee.
- · Ligaments or tendon injury.
- · Quad cramps.
- Knee Sprain.
- Weak knees.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee





AD-414

Patellar Support

Provides enhanced kinesthesia for promoting blood circulation and healing with customized compression and padding.

Indications:

- Patellofemoral pain syndrome.
- Patellar tendonitis and arthritis.
- · Osgood-Schlatter disease.
- Patellar dislocation and subluxation.
- Mild instabilities of the knee joint.

Size	Inches	Cm
Universal	12-20	30-50



*Measure circumference around the knee




Compression Stocking Mid Thigh

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference a mid thigh approx 6" above knee




AD-704

Compression Stocking Below Knee

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

Size	Inches	Cm
S	8.5-11.0	21.5-27.90
М	11.0-13.50	29.90-34.20
L	13.50-16.0	34.20-40.60
XL	16.0-18.50	40.60-46.90
XXL	18.50-21.0	46.90-52.50



*Measure circumference around calf approx 6" below knee





AD-710

Medical Compression Stocking Mid Thigh -Classic

Medical compression stockings are specially designed to apply pressure to your legs, helping to maintain blood flow and reduce discomfort and sewlling.

Indications:

- Maintain blood flow
- Reduce discomfort & swelling
- Improving lymphatic drainage
- Managing orthostatic hypotention

Perimeter/ Model	S	М	L	XL	1
сВ	20-22	23-25	26-29	30-32	AI -
сС	30-39	33-41	36-45	39-48	00
cD	29-38	32-40	35-44	38-47	ec (=====
сС	49-63	54-66	58-71	60-75	· ·

*Measure circumference a mid thigh approx 6" above knee



Medical Compression Stocking Below Knee -Classic

Medical compression stockings are specially designed to apply pressure to your legs, helping to maintain blood flow and reduce discomfort and sewlling.

Indications:

- Maintain blood flow
- Reduce discomfort & swelling
- Improving lymphatic drainage
- Managing orthostatic hypotention

Perimeter/ Model	S	М	L	XL
сВ	20-22	23-25	26-29	30-32
сC	30-39	33-41	36-45	39-48
cD	29-38	32-40	35-44	38-47
сC	49-63	54-66	58-71	60-75

*Measure circumference around calf approx 6" below knee



Supporting Health



AD-406

Knee Immobilizer (14")

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

Indications:

- · Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- · Acute sprains/strains of the knee.
- · Post-cast care.
- Post/Preoperative care.

Size	Inches	Cm	
S	14.8-17.2	37-43	,
М	17.2-19.6	43-49	6
L	19.6 - 22	49-55	
XL	22-24.4	55-61	_
XXL	24.4-26.8	61-67	






Knee Immobilizer (19")

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

Indications:

- Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- · Acute sprains/strains of the knee.
- · Post-cast care.
- · Post/Preoperative care.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
Ĺ	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee




AD-408

Knee Immobilizer (22")

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

- Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- Acute sprains/strains of the knee.
- Post-cast care.
- Post/Preoperative care.

Size	Inches	Cm
S	14.8-17.2	37-43
М	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx $6^{\prime\prime}$ above knee





AD-411

R.O.M Knee Brace (22")

Our R.O.M knee brace provides locking at any angle or motion between any two angles with intervals of 15 degrees.

Indications:

- Arthritis/Osteoarthritis.
- Hamstring Injury.
- · Knee ligament instability.
- Knee sprain.
- Lateral Collateral Swelling & inflammation.
- Total knee replacement.

Size	Inches	Cm
Universal	14.8-22	37-55



*Measure circumference around mid thigh approx 6" above knee




Calf Support

Hypoallergenic elastic fabric from unique fibres blended with cotton prevents skin irritations and perfect compression.

- Weak or injured calves, fatigued legs, and stressed muscles.
- Aching calf muscles.
- Cramps and painful calf conditions.
- Muscle strain.

Size	Inches	Cm
S	9.8-11.3	25-29
М	11.3-12.8	29-33
L	12.8-14.3	33-36
XL	14.3-15.8	36-40



*Measure circumference approx 2" below knee

AD-417

FOOT, HEEL & TOES SUPPORT





AD-413

AD-706

Foot Drop Splint (Right/Left)

Useful for patients under critical care of foot injury. Lightweight, robust and flexible polypropylene shell in one piece.

Indications:

- Muscle weakness due to stroke, spinal cord injury, muscular dystrophy, cerebral paralysis, poliomyelitis, multiple sclerosis and peripheral neuropathy.
- · Recovery from surgical or neurological treatments for foot drop.
- Muscle tightening and runner's cramp.
- Plantar fasciitis.
- Achilles tendonitis.

Arch Support

Size	Europe	USA	ММ
СН	25-30	8-12	157-190
S	34-36	2.75-4.35	217-230
М	37-39	5.25-6.75	237-250
L	40-42	7.50-9.00	257-270
XL	43-45	9.75-11.5	277-291

*Measure Shoe Size



Silicone based material absorbs shock, reduces friction and improves general foot comfort. Controls foot position and prevents arch pain.

Indications:


- Flat foot and associated problems.
- · Walking discomfort.

Toe Separator

and realign the toes.

• Plantar fasciitis.

Size	European	American
S	29-34	CH 11-AD 3
М	34-39	3.0-7.0
L	39-43	7.0-10.50



*Measure Shoe Size




Indications:

· Bunion, corns and calluses pain.

Odour-free, silicone based toe separator helps release pressure

- Hallux Valgus.
- Toe misalignment.
- Plantar Fasciitis.
- Crooked toes, claw toes & hammertoes.

Size	European	American
S	34-36	2.8-4.4
М	37-39	5.3-6.8
L	40-42	7.5-9.0



*Measure Shoe Size





AD-702

Insole (silicone)

Made of medical grade silicone for durability and comfort.

Indications:

- Achilles Tendonitis.
- Heel pain.
- · Heel spurs.
- Pain due to walking.

Size	European	American
S	34-36	2.8-4.4
М	37-39	5.3-6.8
L	40-42	7.5-9.0
XL	43-45	9.8-11.5



*Measure Shoe Size



AD-701


Heel Cushion (silicone)

Made of medical grade silicone for durability and comfort.

Indications:

- · Achilles Tendonitis.
- Heel pain.
- Heel spurs.

Size	European	American
S	34-36	2.8-4.4
М	37-39	5.3-6.8
L	40-42	7.5-9.0
XL	43-45	9.8-11.5



*Measure Shoe Size





AD-707

Metatarsal Pad

Decreasing the pressure on the ball of your foot by supporting the metatarsal bone.

Indications:

- Plantar Fasciitis.
- Plantar Fat Pad Atrophy.
- Metatarsalgia (stone bruise)

Size	European	American
Universal	34-45	2.8-11.75

*Measure Shoe Size



AD-418

Cast Shoe

A lightweight cast shoe which improves balance and helps early ambulation.

- Neuropathic diabetic foot ulcers.
- Forefoot and mid-foot injury.
- Post ligament, soft tissue and tendon surgery.

Size	European	American			
S	34-36	2.8-4.4			
М	37-39	5.3-6.8			
L	40-42	7.5-9.0			
XL	43-45	9.8-11.5			



*Measure Shoe Size



HOT & COLD PACKS





Hot & Colds Pack

Pack made of heating refrigerant gel which instantly relieves muscle pain and other aches due to wounds, sprains, etc.

Indications:

- Lumbago (Low Back Pain).
- Joint pain & swelling.
- Muscle spasms.

Si	ize
Unix	versal

PHYSIO AIDS





AD-1301

Exercise Band (Actiband)

It can be used for a variety of applications. Used particularly lower body exercise to increase strength, body Balance and Strengthen muscles.

- Improve strength & flexibility in glutes, Knees and hips.
- Convenient and Portable.
- Powder-free surface for better gripping.

Resistance Level	Colour Code	Resistance in Lbs*	Resistance in Kgs*
Light	Yellow	09.98	4.5
Medium	Red	11.86	5.4
Strong	Green	13.73	6.2
Extra Strong	Blue	15.60	7.0
Super Strong	Black	16.85	7.7

OTHER PRODUCTS



Size cm

6 cm

8 cm

10 cm

15 cm



Disposable 3 Ply Face Mask

3-ply mask is a disposable mask that has 3 layers of non-woven fabric to help you breathe clean air.

Size	
Universal	

Product Code

AD-1402

AD-1403

AD-1404

AD-1405

Indications:

- Adjustable nose bridge
- · Seamless crimping & breathability
- Soft non-woven earloops



Cotton Crepe Bandage B.P. (Acticrepe)

Cotton Crepe Bandages provide optimum support to the injured area and helps in relieving pain.

٠							
п	n	\sim	10	at	-	n	c.
п	ш	ıu	16	aι	ıv		э.

- Helps relieving muscle strain and joint pain.
- Offers firm support and compression along with warmth.
- Made of skin-friendly material to ensure comfort.
- Material is soft and absorbent.
- · Reusable and washable.



Supporting Health

NOTES





NOTES



Supporting Health



















CONTACT SALES

Phone: +91 141 4015880

Email: sales@adoremedical.com **Website:** www.AdoreMedical.com

Address:

Adore Medical Private Limited Plot No. A-189(A), Road No.01-D, Vishwakarma Industrial Area, Jaipur - 302013, Rajasthan, India









